



Bolling Mitral Valve Retractor

Uncompromised Exposure


BOLLING RETRACTOR

- Better Visibility
- Increased Working Space
- Safe and Stable Retraction



Is your retractor in the way?

Bolling Mitral Valve Retractor




Improve Exposure with the Bolling Mitral Valve Retractor:

- Low Profile Handles are never in the way
- Frame and handles do not interfere with sutures
- Blades stay in place and are easy to set up
- Provides exposure for a variety of procedures:
 - Trans Septal Approach
 - Direct Atrial Approach
 - Mitral Valve Repair
 - Mitral Valve Replacement
 - MAZE



Don't let the competition get in the way.

Bolling Mitral Valve Retractor




- Attaches directly to most spreaders
- Lowest profile Mitral Valve Retractor
- Spreader not included (Morse Spreader 5768)



3 frame options

Bolling Mitral Valve Retractor



Bolling Retractor 91070CV

Qty	Item#	Description	Length	Qty	Item#	Description	Length
(A) 1	43100	Malleable Paddle	5"	1	43110	Medium Basket	6.5"
1	43101	Malleable Paddle	6.5"	1	43111	Medium Basket	8"
1	43102	Malleable Paddle	8"	(E) 1	43112	Large Basket	5"
(B) 1	43103	Cushing Blade	5"	1	43113	Large Basket	6.5"
1	43104	Cushing Blade	6.5"	1	43114	Large Basket	8"
1	43105	Cushing Blade	8"	(F) 2	47200	CV Clamp L-Bar	2"
(C) 1	43106	Small Basket	5"	(G) 1	47201	CV Clamp T-Bar	8"
1	43107	Small Basket	6.5"	(H) 1	47202	CV Clamp L-Bar	6" x 8"
1	43108	Small Basket	8"	(I) 3	47204	Universal Joint	
(D) 1	43109	Medium Basket	5"	1	50000CV	Instrument Case	15" x 9.5" x 1.5"

Thompson 
SURGICAL INSTRUMENTS, INC.

10170 East Cherry Bend Road
Traverse City, MI 49684
231.922.0177 • Fax 231.922.0174
www.thompsonsurgical.com

FREE TRIAL
Put it to the test
800.227.7543



jcb01ss0209